Pharma International

Das unabhängige Informationsorgan für die **Welt-Pharma & Biotech** Industrie
The Independent Informational Organ for the **World's Pharmaceutical & Biotech** Industry

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A Well-Rounded Invention

A new pharmaceuticals feeder
for bulk ingredients with innovative
design

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A Well-Rounded Invention

by Rolf Welsch

A new pharmaceuticals feeder for bulk ingredients with innovative design

If feeding equipment is to be used for pharmaceutical ingredients in the production of drug products, it must comply with the strict regulations of high purity production. Moreover, the applicable regulations and guidelines of the Pharmaceutical Industry must also be complied with, using, of course standardized specific components wherever possible. Since one basic requirement is highest precision and reproducibility in order to guarantee quality, reliable and uniform flow of the bulk material in the feeders must be guaranteed.

Due to many years of experience in the Foods Industry with the development of hygienic design versions of feeding units based on the FlexWall® Principle, it was logical to expand this globally known mass flow feeder to include applications in the pharmaceutical industry. Implementation of constructive information from the USDA and compliance with the 3A Standards had already led in the Hygienic Design to FDA-acceptable parts, which had even been integrated in the standard program.

The following text presents a new patent-pending FlexWall® Feeder which meets the even more stringent conditions in the pharmaceutical industry. This presentation shows the construction of the feeder, the cleaning possibilities and provides information on its areas of use. Finally, there is discussion of certain unique control features.

The basic requirements for the unit were clearly defined at the start of development.

It had to

- be so constructed that it is easy to dismantle and clean.
- be formed inside to prevent any dead areas and corners where residues could collect and cause contamination.
- have a novel outer form with a rounded surface to prevent adherence and accumulation of bulk substances and cleaning fluids.
- be sealed preventing ingredient leaking out and contamination seeping in.
- be guaranteed that only «FDA approved materials» are used.

Figure 1 shows the new design:

The first requirement cited easy to dismantle and clean which meant reducing the connecting parts to a minimum and eliminating the the need for tools to dismantle the feeder. This in turn meant the classical stirring agitator couldn't be used. The state of the art principle of aseptic sealing flat joints was to be applied in any event.

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Well combined

Figure 2 shows the newly developed hopper form of the FlexWall® Feeder. The basis for the invention is clearly shown. The Figure shows in Pos. 2 that the inlet opening of the feeding hopper has a reinforced elastic edge sealant which is affixed to the hopper as an «integrated» gasket. The same solution was created for the horizontal flexible hopper outlet and the connection to the screw mount. Figure 2, Pos. 4 and 7, shows the reinforced elastic edge for these hopper openings.

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Fig. 2: Flex®Wall-Hopper with the «integrated» gaskets

Figure 3, Pos. 1 and 2 shows in detail how the hopper inlet and outlet openings lie loose in the flanged connection.

This edge construction has two functions. On the one hand, it serves to attach the hopper to the housing and on the other, it seals the inside of the feeder, that is the product space, optimally against outside influences. Thus, no separate sealant is needed in addition to the hopper mount. The gasket is an integral component of the flexible hopper, so it is impossible that the gaskets jam or shift during assembly. The feeder extension hopper can be safely positioned anytime with this solution and the slightly modified flanged connection.

These «integrated gaskets» are preferably round and fit to the inner side of the flanged connector. This creates hygiene compliant aseptic sealing. It corresponds to the aseptic clamp tube connectors used in the pharmaceutical industry which are approved for production of drug products and meet the pertinent standards and FDA requirements (Figure 4).

In assembling, the reinforced elastic edge of the hopper body is placed loosely in the corresponding groove and the part to be mounted fitted with its groove. This automatically results in coaxial centering. The clamping ring is stretched so that the front sides of the flanges abut. This produces two effects:

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Fig. 3: Detail of the front of the hopper with inlet and outlet openings

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Fig. 5: Top view of the feeding hopper, there are no disrupting blades

- The elastic edge is reshaped inside to create a smooth, dead space free passage. Penetration of bulk ingredient from inside or of dust particles from outside into the area between the elastic edge of the hopper and the groove is not possible.

- The relatively sharp edge on the inner side of the flanged connector exerts a very strong pressure on the elastic edge of the hopper, guaranteeing very high tightness both with underpressure and overpressure in the hopper.

During operation, the gasket should maintain the area pressure attained on assembly over a wide range of temperatures. During assembly, the PUR-gasket changes shape very little and thus adapts exceptionally well to micro and macro roughness of the gasket surfaces. The swelling value, pressure shape changes and hardness changes are negligible with respect to many usable cleansers.


Abb. 4: Detail der Trogverbindung

Fig. 4: Detail of hopper connection

In the best form

Optimization with respect to avoiding dead space had already been achieved by the optimal inner shape of the FlexWall® Feeders:

- Its interior construction itself has few crevices. Since the paddles for mass flow activation are outside, there is no break through to the gears, as is required for stirrers. There are also no cracks between the sheath of a removable stirrer its gears plus a few corners and edges on the blades. The top view (Figure 5) of the feeder hopper shows clearly the free passage down to the feeding screw (Figure 5, Pos 3).
- There is also no transitional edge to a lower screw zone (Figure 2, Pos 1). These are necessarily the result of construction when stirring agitators are used, because screws can only be placed under the blades.
- The screw connector is placed so far back in the flexible hopper that there is automatic relief from the pressure of the bulk substances (Figure 2, Pos 5). This naturally reduces any microbiological exposure of pharmaceutical bulkware due to too-intensive contact with the mounts on the gearing side.

In principle the Polyurethane surface is largely homogeneous and pore-free, and is comparable to the steel surfaces mentioned later with respect to contamination. The mechanical resistance capacity of PUR is known to be even higher than that of steel. Additional safety in using abrasive bulk substances!

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Grab the opportunity

Due to a completely new development, there was the opportunity to change the outer shape to fully meet the requirements in pharmaceutical production. Both housing halves were made cylindrical or tapered so that no bulk materials can accumulate (Figure 6). Cleaning water is quickly and completely drained. Since the motor in the lower part of the feeders is completely surrounded by the housing, use in a clean room is permissible.

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Fig. 6: Top view of the feeding hopper, there are no disrupting blades

Usual constructions, on the other hand, consist of housing constructed around existing mechanics and have horizontal surfaces and screw connected coverplates with additional gaskets. Even a welded connection of the hopper to the funnel may have disadvantages if the unit is so big that visual inspection is not possible everywhere. The degree of cleanliness has to be determined by the swab test. And the hard-to-reach areas have to be defined in the validation.

The specification of good stainless steel surface quality varies among the pharmaceutical users. In some applications (frequently in the USA), stainless steel 304/304 L (corresponding to DIN 1.4301/1.4306) is considered adequate for parts which have no direct contact with the bulk substances. For reasons of general acceptance in the pharmaceutical industry, we have used stainless steel Type AISI 316 L (DIN 1.4404) even for the housing of the FlexWall®-Feeder described here, although this does not come into contact with the product (Certificate DIN 50049 2.2 or 3.1.B on request). The sheets have a roughness depth of $Ra \leq 0.8\mu\text{m}$ for product contact parts and $Ra 1.2 \leq \mu\text{m}$ on the outside. They are directional-ground to enhance draining of cleaning agents. For product contact parts, a roughness depth of only $\leq 0.4\mu\text{m}$ is also available on order. Aluminum parts have been avoided as far as possible since some cleaners are aggressive to surfaces. The soldered seams are completely soldered, inert and free of foreign matter and also have the required roughness depth goes without saying.

Every gram counts!

From the process technical point of view, thermic and mechanical influences on the sensitive pharmaceutical bulk ingredients, such as occur in stirring, are also avoided by applying the FlexWall® Principle. Expensive additional maintenance which is always required by the stirrer is avoided when using the FlexWall design.

The use of twin screws with diameter of only 12 mm and the smallest gradient meets the requirement for small

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production starting at 50 g/h in tablet pressing and extrusion. The needs of the pharmaceutical industry are adequately covered with further equipment types and spiral screws of up to 40 mm in diameter. One special advantage is that precisely fitting flexible hoppers have been developed to fit the various screw diameters. This was the only way to fulfill the known demand for minimum residual quantities of the often expensive pharmaceutical powders. In such operations, every gram counts!

Since many uses in the pharmaceutical industry require feeding of fine, dusty powders, it made sense to consider the possibility of explosion prevention in compliance with ATEX 95 right from the start. None of the new design characteristics had to be sacrificed to meet these requirements.


Abb. 7: Flexible Verbindung

Fig. 7: Flexible sleeves

Provision is made for cooling the motor and providing an protective gas connection for the dynamic sealing (both optional). The operator will use laboratory pressurized air of ISO 8573.1/class 1 quality whenever possible for pharmaceutical applications. If this is not the case, it is recommended to produce this using special oil-free hygienic laboratory compressors. These have high noise reduction (61 dB), usually also contain a pressure reducer, microfine or active charcoal filters and sometimes dryers. The quantity of air needed to cool the feeder motor is only 60 l/min., so low that the smallest available laboratory compressor is sufficient.

Flexible sleeves should be used to connect the Pharma Feeders as loss in weight feeder to filling equipment or conveyors which are provided with aseptic clamp rings like the hopper. They are offered with the unit (Figure 7). It is frequently recommended to have additional air refeed pipe from the discharge tube to the lid of the weighing container, since ventilation to the outside is then not necessary. The number of parts to be cleaned is then greater, but it saves problematical ventilation filters in the clean rooms (Figure 8).

Clean and tidy

The stream of the spray ball hits against the flexible high resistance polyurethane hopper of a FlexWall® Feeders,

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Fig. 8: Feeder with refeed pipe

making it vibrate (hydromechanical factor). This produces a better cleaning effect than on steel surfaces. The installation of such a spray ball is seen in Figure 9. Many users then decide against using rotating spray balls, so that they can be completely constructed of stainless steel. For a

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Fig. 9: Installation of a cleaning jet

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FlexWall® Feeder, the construction with integrated spray balls is also sufficient, since there are no spraying shadows caused by any blades.

Only slight soiling of the dynamic sealing at the gearing block is to be expected when the optional barrier gas connection is in place. However, it can be easily removed for cleaning by releasing one clamp ring. Most parts of the Pharma-Feeders described here can be placed without problems in autoclaves for sterilization, taking temperatures and exposure time into account. But remember that this cleaning has to be clearly described for validation.

...., Control is better

The user might think that the use of volumetric feeders would simplify cleaning, since additional parts of a weighing technique are omitted. But remember that continuous weighing saves frequent calibration on the open bulk flow (dust!) and there is no need for an additional filling level probe, which also make cleaning a volumetric feeder more complicated.

An in-process weight control in gravimetric feeding is also no problem in a clean room. Control of the preweighed mass must, however, be possible without the use of sample vessels. For this, the balance control has a sampling button. When activated, it measures the continuous feeding flow over time of, for example, 1 minute during production (online) and compares the theoretical value to the actual value. Usually, the user evaluates the results according to NAMUR (statistical 2-Sigma-Method).

Reliance on the weighing technical precision of the loss in weight feeder can justifiably be very high, since weight control is possible at any time by using a calibration weighting stone on the scale, which is located in the cleanroom/isolator. The high precision of a gravimetric system leads to a reduction in the quantity of expensive pharmaceutical bulk ingredients used, thanks to the tighter theoretical value tolerances.

In using feeders, there are unfortunately a few working steps which cannot be rationally automatized. Assembly in an isolator, for example, must be done through the gloves. This applies to precleaning of parts, changing the screw, etc.

Isolators often have an integrated manual vacuum hose with sprayer for intensive flushing of the unit. The light flexible FlexWall® Hopper can be optimally cleaned with this by vacuuming off the slight residuals then flushing. The hopper can be detached in the isolator at any time without damaging other parts preventing it from forming grooves, as may happen in using metal hoppers.

Even when working through the gloves, cleaning is very easy. Doing away with the agitator, as mentioned several times, simplifies assembly and dismantling to a great degree. Sharp edged blades, on the other hand, can easily destroy the gloves!

Validation

Connection to the customer's Controller, who undertakes administration and storing of the measured data, is made

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via interfaces on a balance controller (BUS, Ethernet etc.). The data from the feeder controller are thoroughly processed to be adequate or as required for the customer's validatable software (e.g. in compliance with the assessment criteria of FDA CFR 21 Part 11). For the pharmaceutical customer, the feeder controller is a sort of Black Box, which is specified as follows (abbreviations are in Appendix).

- The feeding performance and required precision are checked in a **PQ**.
- The installation site, type of filling, EX-zone divisions and type of cleaning are determined in the **FS**.
- The description of the **HS** means for the balance supplier that he must confirm compliance with the regulations (e.g. FDA /USDA, Customer's specifications) as well as appropriate design and construction. The documents necessary for the accepted parts, such as clamps etc. are attached.
- **OQ** can mean for the balance that the function of the most important working points is checked on the basis of the measuring protocols.

The password access required by the FDA is, of course integrated in the balance controls, and the calibration required is possible and described precisely. The user has the option of installing only the operation control Type OP 1 with very small dimensions (144x144x30mm, HxWxD) in the processor machine. Feeding control is often integrated in the processing machine these days, but for users in the pharmaceutical industry, assembly of the motor and regulators outside the cleanroom is better to avoid needing expensive pharmaceutical cabinets for these parts.

PQ:	Performance Qualification	FDA:	Food and Drug Administration
FS:	Functional Specifications	URS:	User Requirements Specification
HS:	Hardware Specification	GAMP:	Good Automated Manufacturing Practice
OQ:	Operation Qualification	ANSI:	American National Standards Institute
CGMP:	Current Good Manufacturing Practice	EHEDG:	European Hygienic Engineering & Design Group
FAT:	Factory Acceptance Test	URS:	User Requirement Specification

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Summary

The new FlexWall® Feeder in «Pharma Design» presented here has already in use by various pharmaceutical producers. It is important for the pharmaceutical customer to find a supplier of feeding equipment who can understand not only the handling of bulk ingredients, but also the problems which arise in complying with all the regulations in addition to the customer's own specific requirements.

The supplier should be ISO-certified to be able to produce the Quality Plans in support of his own validation process. The extensive documents should also be offered, at least as an option. Pharmaceutical machines are safely packed and stored separately in certified construction to prevent any later damage to surfaces. ■

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